

Brief instructions



1 Remove the foil cover from the cuvette. The liquid surface should lie between the two lines marked on the cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette.



2 Fill the the orange striped capillary with the sample (10 µl) up to the white stopper. Ensure that there are no air bubbles in the capillary. Wipe away any residual sample from the outside of the capillary.



3 Place the sample (10 µl) into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.



4 Close the cuvette tightly with wrCRP reagent cap. Do not press down the turquoise inner part of the cap. Once the sample has been added to the buffer the assay should be run within two hours.



5 Choose Measure on the display of the QuikRead go instrument.



6 Insert the cuvette in to the measurement well of the instrument. The bar code should be facing you. The display shows how the measurement is proceeding.



7 When the measurement is completed, the result is shown on the display and the cuvette will rise up from the measurement well automatically. With the QuikRead go wrCRP+Hb test the Hb result is shown below the CRP result.

Please consult the instructions for use before performing the test.

Reagent storage



Kit box

Storage of unopened kit at 2–25°C: until the expiry date of the kit



Unopened cuvettes after opening the foil pouch

Storage at 2–8°C: 6 months

Storage at 18–25°C: 3 months



Reagent cap tube

Unopened tube at 2–25°C: until the expiry date of the kit

After the first opening at 2–25°C: 6 months



Opened cuvette

The sample is stable in the buffer for two hours. Perform the test within **two hours** from opening the cuvette. Let the buffer warm up to room temperature before performing the test.



wrCRP Controls (see table below)

Unopened vials at 2–8°C: until the expiry date

Opened vial at 2–8°C: up to two months



QuikRead go Hb 10 µl Control (see table below)

Unopened vial at 2–8°C: until the expiry date

Opened vial at 2–8°C: up to one month



Orion Diagnostica Oy
Koivu-Mankkaan tie 6 B, P.O.Box 83, FI-02101 Espoo, Finland
www.oriondiagnostica.com, www.quikread.com, www.tacklearm.com

8046-05EN, 11/2018. International version. Not intended for the US market.

Product	Cat. no.
QuikRead go wrCRP, 50 tests » with capillaries (10 µl) and plungers	146521
QuikRead go wrCRP, 500 tests	146522
QuikRead go wrCRP+Hb, 50 tests » with capillaries (10 µl) and plungers	146255
QuikRead go wrCRP Control 1 ml	153038
QuikRead go wrCRP Control High 1 ml	153039
Ultrasensitive CRP Control 1 ml	68257
QuikRead go Hb 10 µl Control 1 ml	153656
QuikRead go 10 µl capillaries, 50 pcs	147851
QuikRead Plungers, 50 pcs	67966
QuikRead go Instrument	133893

Find more information: www.quikread.com